ANI Monitor V2

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STATE of the ART

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	Design of the study	n	Results	Conclusion	References
Awake patients	During labor ANI and NRS during and between uterine contractions	45	Negative linear relationship between ANI and NRS ($p=-0.18$)	Good correlation between ANI and NRS during labor	Le Guen M et al., Int J Obstet Anesth, 2012.
	ANI and NRS in PACU after general anesthesia, ANI performance for the pain detection (NRS > 3) in PACU	200	Negative linear relationship between ANI and NRS ($r^2=0.41$) Good pain detection performance in PACU (AUC ROC=0.86) ANI ≤ 57 Se =78% and Sp=80%	Good correlation between ANI and NRS during immediate postoperative pain	Boselli E et al., Br J Anaesth, 2013.
	ANI variations and VAS scores during physical therapy procedure, 24 and 48 hours	12	Decrease of ANI when VAS > 30. Good pain detection (AUC ROC=0.76, ANI Se=76%, SP=78%)	ANI can be used for pain management during physical therapy	De jonckheere J et al., Conf Proc IEEE Eng Med Biol Soc, 2014.
	ANI variations during the projection of a violent film, comparisons between controls and anorexia nervosa patients	24	Decrease of ANI during the projection, with a delay anorexia nervosa patients	ANI is decreased after an unpleasant emotion	Rommel et al. Psychiatry Res, 2015.

	Design of the study	n	Results	Conclusion	References
Anesthetized patients	Pediatric (8 ± 5 years old) ANI variations after tetanic stimulation (50 mA, 50 ms, 5s) with desflurane and different concentrations of remifentanyl (0.04 to 0.2 µg/kg/min)	12	Decrease of ANI after tetanic stimulation, more important with the lowest remifentanyl concentration	ANI is decreased after tetanic stimulation during pediatric surgery and seems more sensitive than the skin conductance	Sabourdin N et al., Paediatr Anaesth, 2013.
	Pediatric (2 – 16 years old) Locoregional anesthesia before surgery and sevoflurane only, ANI performance for the detection of a locoregional anesthesia failure (increase HR ≥ 10 % 2 min after incision)	58	ANI performance for the detection of locoregional anesthesia failure: AUC ROC=0.75 (0.61 - 0.88) ANI ≤ 51 , Se=79% and Sp=62%	ANI is able to detect locoregional anesthesia failure	Migeon A et al., Paediatr Anaesth, 2013.
	Intravenous anesthesia with propofol and remifentanyl ANI variations after tetanic stimulation (50 mA, 60 Hz, 30 s) at different targets of remifentanyl (0., 2 and 4 ng/mL)	25	Decrease of ANI after tetanic stimulation	ANI is efficient for the detection of pain response during surgery under general anesthesia	Gruenewald M et al., Br J Anaesth, 2013
	One hundred and eighty adults presenting for elective surgery were included. All received total intravenous anaesthesia with propofol adjusted to entropy and remifentanyl adjusted to the ANI.	180	Anaesthesia was achieved without any episode of reactivity in 160 (89%) patients. The median remifentanyl dose was 0.042 [0.040-0.044] µg.kg ⁻¹ .min ⁻¹ . At 24 hours, the maximal NRS pain score was 2 [2,3]. One hundred and fifty-five patients (86%) did not receive any postoperative opioids.	ANI can be used to adequately guide intraoperative remifentanyl administration during vascular surgery. Such guidance resulted in low remifentanyl consumption, low postoperative pain rates and low opioid rescue analgesia.	G Daccache et al. Anaesth Crit Care Pain Med. 2016

	Design of the study	n	Results	Conclusion	References
Anesthetized patients	General anesthesia under sevoflurane and remifentanyl ANI variations after tetanic stimulation (50 mA, 60 Hz, 30 s) at different targets of remifentanyl (0, 2 and 4 ng/mL),	24	delta ANI significantly indicated patient's movement after tetanic stimulation (Se=77% and Sp=84%)	ANI is efficient for the detection of pain response during surgery under general anesthesia	Gruenewald M et al., Minerva Anesthesiol, 2014.
	Intravenous anesthesia with propofol and remifentanyl Abdominal laparoscopic surgery ANI variations at different times and after tetanic stimulation (80 mA, 100 Hz, 5 s)	15	Decrease of ANI during different nociceptive stimuli.	ANI is able to detect painful stimulations under general anesthesia	Jeanne M et al., J Comput, 2012.
	Intravenous anesthesia with propofol and sufentanil bolus Knee arthroplasty ANI performance for the detection of perioperative hemodynamic reactivity (increase > 20% HR and/or BP in 5 min)	27	Good performance for the hemodynamic reactivity (AUC ROC=0.92) ANI ≤ 63 Se=80% and SP=88%	A drop of ANI values predicts a hemodynamic response lead by pain	Jeanne M et al., Clin J Pain, 2014.
	TIVA propofol and remifentanyl Laryngoscopy suspension procedures ANI performance for the prediction of hemodynamic reactivity and sedation.	50	Good performance for the hemodynamic reactivity (AUC ROC=0.88) ANI ≤ 55, Se=88% and Sp=83%	The prediction of hemodynamic reactivity is slightly reduced by the use of sevoflurane and/or fentanyl	Boselli et al. Minerva Anesthesiol, 2015.
	General anesthesia with sevoflurane and fentanyl General or orthopedic surgery ANI variations at different times PK for the prediction of increase HR and Systolic BP > 10%	30	Decrease of ANI after nociceptive stimulation Increase of ANI after fentanyl administration Modest prediction probability for increase HR (PK=0.61), increase BP (PK=0.59)	The prediction of hemodynamic reactivity is slightly reduced by the use of sevoflurane and/or fentanyl	Ledowski T et al., Acta Anaesthesiol Scand, 2013.
	General anesthesia with halogen and remifentanyl ORL surgery or inferior limb orthopedic surgery ANI performance before extubation for the prediction of immediate postoperative pain (NRS >3)	200	Negative linear relationship ANI before extubation and NRS in PACU (r ² =0.33) Good ANI performance for the prediction of NRS > 3 (AUC ROC=0.89) ANI ≤ 50, Se=86%, Sp=86%	ANI values at the end of the surgery are able to predict postoperative pain	Boselli E et al., Paris:SFAR, 2014.
	General anesthesia with sevoflurane and fentanyl Abdominal hysterectomy procedures Effect of 0.5 mg/kg ketamine administration on ANI	20	No modifications of ANI 5 min after ketamine administration	Ketamine administration does not interfere with ANI response during surgery	Bollag L et al., J Clin Monit Comput, 2014.
	This study compares the analgesic indices Analgesia Nociception Index (heart rate variability), Surgical Pleth Index (photoplethysmography), and pupillary dilatation, to heart rate, mean arterial pressure, and bispectral index, with regard to diagnostic accuracy and prediction probability for nociceptive response.	38	Under propofol sedation, sensitivity and specificity of the Analgesia Nociception Index (PK = 0.98) for detecting painful stimulations were high compared to heart rate (PK = 0.74), mean arterial pressure (PK = 0.75), and bispectral index (PK = 0.55)	The Analgesia Nociception Index, the Surgical Pleth Index, and pupillary dilatation are superior in detecting painful stimulations compared to heart rate and mean arterial pressure	Funcke S et al, Anesthesiology 2017

ANI Monitor V2

Spécifications techniques

Général

Parameter	Specification
Power Requirements	100-240 VAC through AC power adapter
Mains Frequency	50/60 Hz
Battery Type	No battery
DC Input	12V 3,4A 40W
Patient Leakage Current	<5 μ A @ 220V and 50 Hz

Connecteurs

Parameter	Specification
DC Input (monitor)	3-pin power connector
ECG IN (monitor)	6-pin female connector for sensor cable connection
Export (monitor)	Sub-D9 connector to export data in real time
Data Export (Monitor)	USB connector to export data and snapshot to USB stick
Sensor cable	6-pin male connector
Sensor (Acquisition Device)	5-pin female connector for sensor

Eport

Parameter	Specification
Export Protocol	UART interface
Data Export	USB interface

Environment

Parameter	Specification
Cooling Method	Convection. Fan less
Temperature Operating Storage	5°C to 40°C -20°C to 60°C
Dimensions Monitor (width x height x depth)	256 x 214 x 81 mm
Weight Monitor	2,5 Kg
Finish Monitor	Front : black and orange Back : white

Écran

Parameter	Specification
Type	Color Liquid Crystal
Size	200 mm (8 inches)
Resolution	800 x 600 pixels
Active Viewing Area	173 x 130 mm
Pixel pitch	0.216 x 0.217 mm

Mdoloris Medical Systems SAS
Biocentre Fleming
Bâtiment C Epi de Soil
270 rue Salvador Allende
59120 LOOS - FRANCE

contact@mdoloris.com

www.mdoloris.com

ANI Monitor V2 is a class IIa medical device, manufactured by Mdoloris Medical Systems,
CE evaluation was performed by Bureau Veritas Italy (1370) . The ANI Monitor V2 name
and logo are registered trademarks. © 2017 Mdoloris Medical Systems. All rights reserved.
MDQUAEN74.ANIV2 v.02