


SUPPLIES & ACCESSORIES

TEMPERATURE



Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
------------	-------------	-------------	-----

Continuous Temperature

Cables

For use with Command Modules 90496/91496, DM3 Monitor 91330 and mCare Monitor 91220



700-0031-00	TruLink Dual Temperature Cable	30 cm/12 in.	1 Each
700-4000-00	Universal temperature probe extension cable, YSI 400 series compatible *Required for use of single use temperature probes/sensors	305 cm/10 ft	1 Each

Reusable Probes/Sensors

Skin Surface



690-0028-01	YSI 700 series compatible with 6.5 mm DIA phone jack type connections	305 cm/10 ft	1 Each
-------------	---	--------------	--------


Esophageal / Rectal









690-0272-00	12 French, YSI 400 series compatible	305 cm/10 ft	1 Each
118-0256-01	12 French, YSI 700 series compatible	305 cm/10 ft	1 Each

Temperature Supplies

A variety of supplies are available through Spacelabs for 90496/91496 modules configured for temperature. A selection of skin, tympanic, esophageal and rectal temperature probes are available in different diameters and lengths.



Order Code	Description	Length/Size	Qty
Single use Probes/Sensors - Requires YSI 400 series cable (700-4000-00)			
Skin			
	701-4000-00 Skin Temperature Sensor		20/Box
Esophageal/ Stethoscope			
	701-4009-00 9 French, YSI 400 series compatible		20/Box
	701-4012-00 12 French, YSI 400 series compatible		20/Box
	701-4018-00 18 French, YSI 400 series compatible		20/Box
Esophageal/ Rectal			
	701-4109-00 9 French, YSI 400 series compatible		20/Box
	701-4112-00 12 French, YSI 400 series compatible		20/Box
Foley catheter with temperature probe			
	701-4214-00 14 French, YSI 400 series compatible		20/Box
	701-4216-00 16 French, YSI 400 series compatible		20/Box
	701-4218-00 18 French, YSI 400 series compatible		20/Box
Tympanic temperature sensor			
	701-4300-00 Adult, YSI 400 series compatible		20/Box
Exergen			
<i>For use with Qube Mini (91389), Qube (91390), Xprezzon(91393) & DM4 (91331) Monitors</i>			
	010-2157-00 Temporal Artery Thermometer		1 Each